

### Cookies

This site uses cookies to offer you a better browsing experience. Find out more on [how we use cookies and how you can change your settings](#).

I accept cookies

I refuse cookies



Live, work, travel in the EU

## COVID-19 In Vitro Diagnostic Devices and Test Methods Database

[Home](#) > [COVID-19 In Vitro Diagnostic Medical Devices](#) > [COVID-19 In Vitro Diagnostic Medical Device - detail](#)

### COVID-19 In Vitro Diagnostic Medical Device - detail

#### Exdia COVID-19 Ag

##### Manufactured by Precision Biosensor, Inc.

Manufacturer website	<a href="http://www.precision-bio.com/en/">http://www.precision-bio.com/en/</a>
CE Marking	Yes
Detection Principle	ImmunoAssay-Antigen
Format	Rapid diagnostic test
Target	Antigen
Commercial Status	Commercialised
Last Update	05/11/2020
<hr/>	
Additional Information	<a href="http://precision-bio.com/en/news/precision-biosensors-covid-19-rapid-Antigen-test-received-ce-mark">http://precision-bio.com/en/news/precision-biosensors-covid-19-rapid-Antigen-test-received-ce-mark</a>
Test Type Result	qualitative
Specimen	nasopharyngeal swabs
Antigen IgG	no
Antigen IgM	no
Antigen IgA	no
Time (Min)	20
Sensitivity Percent	93.9
Specificity Percent	98

The database contains publicly available In Vitro Diagnostic Medical Devices for COVID-19 and it is being updated periodically. Please note that additional performance (as retrieved from manufacturers web pages) is provided only for devices commercially available with CE-IVD mark.

[Acknowledgements](#)

[COVID-19 Test Methods and Devices](#)

[Contact us](#)

[Acknowledgements](#)

This site is managed by the Joint Research Centre

More information on:

[Live, work, travel in the EU](#)

**European Commission**

[Contact the European Commission](#)

[Language policy](#)

[Follow the European Commission on social media](#)

[Cookies](#)

[Resources for partners](#)

[Privacy policy](#)

[Legal notice](#)